

Headspace Vials

CRONUS® 20ml Headspace Vial

Clear Glass, Round Bottom, 22.5x75mm
Part no: V2H20CR-100, Qty: 100 Pack

This vial works with any headspace autosampler. It features a bevelled lip to ensure a reliable seal with septa or stoppers. Round bottom vials are great for samplers that move the vial and need to locate them in a tray or oven.

CRONUS® 10ml Headspace Vial

Clear Glass, Round Bottom, 22.5x75mm
Part no: 405 034, Qty: 100 Pack

This is a 22.5x75mm vial with a rounded bottom specifically designed for accurate headspace analysis. It features a bevelled lip to ensure a reliable seal with septa or stoppers.



Caps and Septa

Pre-assembled

CRONUS® 20mm Magnetic Headspace Seal with SiI/PTFE Septa Small Aperture

Part no: VCA-2004MS-100, Qty: 100 Pack



CRONUS® 20mm Headspace Crimp Cap with SiI/PTFE Septa

Part no: VCA-2004CS-100, Qty: 100 Pack



CRONUS® 20mm Bimetallic Magnetic Cap with PTFE/SiI Septa

Part no: VCA-2004BM-100, Qty: 100 Pack



Separate

CRONUS® 20mm Headspace Crimp Cap

Part no: VCC-2001CS-1000, Qty: 1000 Pack



CRONUS® 20mm Septa SiI/PTFE

Part no: VCS-2004-1000, Qty: 1000 Pack



Crimping Tools

Vial Crimper

For use with 20mm crimp caps
Part no: 200 972

Painted, plated and coated for maximum corrosion resistance. Adjustable stop for optimal sealing performance. Textured handle for excellent grip.



Vial Decrimper

For use with 20mm crimp caps
Part no: 500 446

Decrimpers work in a similar way to crimpers, offering a controlled, low effort method for removing caps. This reduces both sample spillage and breakages.

